

DISINFECTANT *Solutions*



Triple Quick™ Disinfecting Cleaner

- Broad-spectrum power kills a wide range of pathogens
- Versatile formula cleans and disinfects for one-step use
- Ready-to-use formula requires no diluting or mixing
- 3 minute or less dwell time for most claims



State®
Since 1911

CARE FOR WORK ENVIRONMENTS®

Triple Quick™

Disinfecting Cleaner

FAST-ACTING SOLUTION

Triple Quick is a ready to use disinfecting cleaner that prevents cross contamination on facility surfaces. It contains biodegradable surfactants, wetting agents and powerful sequestering agents to provide maximum cleaning power. OSHA Final Rule 29CFR 1910.1030 mandates a variety of regulations regarding employee exposure to potentially infectious materials. Triple Quick is considered to be an appropriate disinfectant per the OSHA Bloodborne Pathogens law.

MARKETS

Triple Quick is the perfect disinfectant in a variety of markets. Utilize Triple Quick in hospitals, healthcare facilities, animal care facilities, athletic facilities, and more.

KILL CLAIMS

Triple Quick is highly effective against bacteria, viruses, fungi and mildew including:

- SARS-CoV-2
- HIV-1 (AIDS Virus)
- Staphylococcus aureus
- Pseudomonas aeruginosa
- Poliovirus Type 1
- Rabies Virus
- Many more
- Mycobacterium tuberculosis (BCg)
- Canine Parvovirus
- Salmonella (choleraesuis) enterica
- Hepatitis A Virus (HAV)
- Avian Influenza A strain H3N2 and H9N2
- Athlete's foot

On EPA List N for use against SARS-CoV-2, the cause of COVID-19. Triple Quick has a 60-second dwell time on List N.

COMMON APPLICATIONS

Triple Quick is effective on all hard, non-porous, inanimate surfaces such as:

- Walls
- Tables
- Telephones
- Light switches
- Hampers
- Bed springs
- Sink tops
- Chairs
- Door knobs
- Shower stalls
- Cabinets

DIRECTIONS FOR USE

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier otherwise enter normally sterile areas of the body. For plastic and painted surfaces, spot test on an inconspicuous area before use.

It is a violation of Federal Law to use this product in a manner inconsistent with its label. Before using, refer to the product label for all applications, procedures and directions.

TECHNICAL DATA:

REGISTRATIONS: EPA Registration Number 1839-83-70799

APPEARANCE: clear, pale liquid

FRAGRANCE: Fresh N Clean, Lavender Meadow

WEIGHT PER GALLON: 8.4 lbs/gal

DISINFECTANT TYPE: quaternary ammonium compound

pH: 11.7

SHELF STABILITY: two years

AUTHORIZATIONS: Meets the requirements for use in Federally inspected meat and poultry plants.

Please review the label and SDS on stateindustrial.com for all product directions, precautions and first aid information.

OTHER STATE DISINFECTANT PRODUCTS




Quick Defense™
Ready-to-Use
Disinfectant



NDC™
Neutral Disinfectant
Cleaner



Discover™
Aerosol Disinfectant
and Deodorizer



Quick Defense™
AQ Wipes
Disinfectant Wipes

PACKAGING

Fresh and Clean

125514 1 QT Bottle/CS12
128341 1 GL Bottle/CS4

Lavender Meadow

125513 1 QT Bottle/CS12
128339 1 GL Bottle/CS4



State Industrial Products
5915 Landerbrook Drive
Mayfield Heights, OH 44124
To Order Call: 1-866-747-2229
www.stateindustrial.com

Puerto Rico
Royal Industrial Park, Bldg "M"
Local #5, Carr 869, km 1.5 Palmas
Cataño, P.R. 00962
To Order Call: 787-275-3185

